

# Model L.TEX 8200 Series High Pressure Pump



The 8200 series has been developed to be used not only for the conventional HPLC pumps but also packing the columns. The 8200 series can be also used for delivering super critical fluid by adding the optional cooling jacket to the pump heads for the SFE/SFC. The 8200 series offers wide range of flow rate using the micro step drive mechanism, compressibility correction function with the pressure sensor, time program function, upper and lower pressure limiter function.

L.TEX 8200 series High Pressure Plunger Pump

Model	Flow rate(ml/min)	Max. pressure	Operation
8210	0.01 ~ 20.0	40MPa	100V~240V
8222	0.02 ~ 50.0	30MPa	100V~240V

Tubing size for IN           Teflon  $\phi 3 \times \phi 2$

SUS  $\phi 1/16$

SUS  $\phi 1/8$

Tubing size for OUT       SUS  $\phi 1/16$

SUS  $\phi 1/8$

\*Specifications and appearance are subject to change without any prior notice



2nd Silk Road Heights, 1721 Katakura-cho, Hachioji-shi, Tokyo 192-0914 Japan

URL: <http://l-tex.net/english/>

Tel: +81-90-4452-7824 Email: [sun\\_inter@rc4.so-net.ne.jp](mailto:sun_inter@rc4.so-net.ne.jp)

20120606-1