

Model L.Tex 8301 Micro Feeder



The model L.TEX8301 can use the gas tight syringes of either 250uL, 500uL or 1mL which are available from the market. Different flow rate can be set depending on the type of syringe. A special joint is supplied which is connected to tubing of 1/16" OD or glass capillary as an exit side tubing.

Model L.TEX 8301 Micro Feeder Specifications

Delivery system	Syringe with a ball screw
Syringe volume	Commercial available gas tight syringe with volume of 250µl, 500µl, or 1ml can be set (no syringe is to be supplied with the system.)
Setting flow rate range	1µl ~ 100µl/min
Max. Pressure	Depending on the gas tight syringe (~5MPa)
Dimensions and weight	120(W) x 125(H) x 230(D)mm Approx. 2kg

*Specifications are subject to change without any prior notice



2nd Silk Road Heights, 1721 Katakura-cho, Hachioji-shi, Tokyo 192-0914 Japan

URL: <http://l-tex.net/english/>

Tel: +81-90-4452-7824 Email: sun_inter@rc4.so-net.ne.jp

20120606-1