

Model L.Tex 8350 High Pressure Micro Feeder



Pressure resistance of the gas tight syringes are not enough for some applications. Therefore the model L.TEX 8350 has been designed to be applicable to the pressure up to max. 20MPa with the same accuracy as that of L.TEX8301 by employing the stainless pump head with 500 μ l volume.

Model L.TEX 8350 High Pressure Micro Feeder Specifications

Delivery system	Syringe with a ball screw
Syringe volume	500 μ l
Setting flow rate range	1 μ l~100 μ l/min
Max. Pressure	20MPa
Dimensions and weight	135(W) x 125(H) x 230(D)mm Approx. 2.5kg

*Specifications and appearance are subject to change without any prior notice



2nd Silk Road Heights, 1721 Katakura-cho, Hachioji-shi, Tokyo 192-0914 Japan

URL: <http://l-tex.net/english/>

Tel: +81-90-4452-7824 Email: sun_inter@rc4.so-net.ne.jp

20120606-1