

Model L.Tex 8551 Syringe Pump



The L.TEX8551 has been developed to inject large amount of the sample precisely for a short time in the line of the supercritical fluid system. The L.TEX8551 is controlled by the PC.

Model L.TEX 8551 Syringe Pump Specifications

Delivery system	Syringe with a ball screw
Setting flow rate range	5.0 μ l~300.0ml/min
Max. Pressure	25MPa
Volume of cylinder	100mL
Operation	PC controll

*Specifications and appearance are subject to change without any prior notice.



2nd Silk Road Heights, 1721 Katakura-cho, Hachioji-shi, Tokyo 192-0914 Japan

URL: <http://l-tex.net/english/>

Tel: +81-90-4452-7824 Email: sun_inter@rc4.so-net.ne.jp

20120606-1